

ICT+

isolator
cleaning tool

MAXIMUM CARE FOR MINI ENVIRONMENTS



 **SATURIX**[®]
by FG CLEAN WIPES

A HIGHER LEVEL OF REGULATORY COMPLIANCE

Conventional isolator cleaning tools require extensive handling to attach and remove pads, which takes time and poses a significant risk for cross contamination. Saturix ICT+, with patent-pending designs, not only eliminates the need to change pads by hand, but also provides a complete system for used-pad management with superior overall usability.

It will help any compounding pharmacy be more compliant with USP and ANNEX 1 standards to ensure the clean and sterile preparation of drugs for the safety of patients and employees.

TOUCH-FREE PAD APPLICATION AND REMOVAL

- Clean pads click on the frame, dirty pads click off the frame—all touch-free
- Reduces cross contamination risk while increasing efficiency

ERGONOMIC, EASY-TO-USE HYBRID DESIGN

- Asymmetrical frame maximizes coverage and effectively cleans tight spaces
- Offset yoke on the frame extends reach—easy to clean high areas

99.9% PARTICULATE RESIDUE REMOVAL

- Microfiber cleaning pads remove 99.9% of particles down to .1 micron*
- Mechanically removes disinfectant residues and other harmful particulates from the surface



Touch-free pad attachment and removal



Precision dosing of pads in package



Used-pad storage within the hood

HYBRID FRAME DESIGN
Easy to clean high and tight spaces

OPTIMIZED SURFACE COVERAGE WITH MINIMIZED CARBON FOOTPRINT

Each pad covers 138 ft² (12.8 m²)

- 100% recyclable pad backer further reduces waste (30% of total replacement pad weight)

Removes **99.9%** of Particulate Residue Down to .1 micron*

* Based on 3rd party lab testing

SATURIX ICT+ ISOLATOR CLEANING TOOL

ICT+ Hardware

SKU#	Product Description	Size	Qty/Pack	Packs/Case
IFHTA	Telescopic Handle and Frame Kit, Autoclavable	Handle: 19.3" – 27.5" (49cm – 70cm) Frame: 8.2" x 2.6" (208mm x 66mm)	1	1
IFH19	Aluminum Handle & Stainless Steel Frame Kit, Autoclavable	Handle: 19" (480 mm) Frame: 8.2" x 2.6" (208 mm x 66 mm)	1	1

ICT+ Microfiber Pads

IPDP8S	Validated Sterile Microfiber Pad in Package	8.2" x 2.6" (208 mm x 66 mm)	4	15
IPDB8S	Validated Sterile Microfiber Pads	8.2" x 2.6" (208 mm x 66 mm)	2	30
IPDP8	Non-Sterile Microfiber Pads in Package	8.2" x 2.6" (208 mm x 66 mm)	4	15
IPDB8	Non-Sterile Microfiber Pads	8.2" x 2.6" (208 mm x 66 mm)	4	15
IPDP8S1	Validated Sterile Microfiber Pad in Single Package	8.2" x 2.6" (208 mm x 66 mm)	1	60
IPA30S	Validated Sterile 70% Isopropyl Alcohol	30 ml, 1 FL OZ	4	15



ICT+ isolator
cleaning tool

 **SATURIX**[®]
by FG CLEAN WIPES

2255 Westover Road, Chicopee, MA 01022
+1 (413) 598-8344 • saturix.com

A Filtration Group Company

ICT-001.2